

Dissolution of MCC in NaOH and Urea

Introduction

This protocol dissolves microcrystalline cellulose (MCC) in NaOH and Urea and is followed by coagulation with 2 M acetic acid.

Reagents

- 🔗 10 g MCC
- 🔗 4 g of 0.5 M urea solution
- 🔗 0.6 g of 0.5 M NaOH solution
- 🔗 2 M acetic acid

Equipment

- 🔗 125 mL Erlenmeyer flask
- 🔗 Magnetic stir plate
- 🔗 Magnetic stir bar
- 🔗 Centrifuge
- 🔗 Glass plate

Procedure

1. In a -5°C cold room, mix 10 g of MCC, 4 g of a 0.5 M urea solution, and 0.6 g of a 0.5 M NaOH solution in an Erlenmeyer flask. Stir the solution for 5 minutes to create a slurry.
2. Hold the solution at -5°C for about 5 hours, then thaw with intermittent mixing at 0°C - 4°C to make a transparent solution.
3. Remove insoluble parts of cellulose by centrifuging at 8000 rpm for 30 minutes.
4. Cast the solution onto a glass plate, taking care to spread it out to about 0.2-0.3 mm, and coagulate it immediately in 2 M acetic acid for 5 minutes.