

Introduction

Protocol for the preparation of electro-competent cells for the use of electroporation, a widely used method for transformation of DNA extracts into cells based on electric pulses.

Materials

- E. Coli Strain with Antibiotic resistance
- LB media
- LB + Antibiotic agar plates
- Antibiotic
- Sterile ice-cold distilled water
- Sterile ice-cold distilled 10% glycerol
- 50 mL falcon tube
- Erlenmeyer culture flask
- Ice
- Eppendorf tubes of 1.5 mL
- Shaking incubator
- Bunsen burner
- Spectrophotometer
- Centrifuge (Eppendorf)
- Centrifuge (50mL Falcon)

Procedure

1. Streak *E. coli* strain onto an LB + Antibiotic agar plate to obtain single colonies and incubate at 37°C **overnight**.
2. Inoculate 2 mL of LB medium + Antibiotic with a single colony of freshly grown *E. Coli* and incubate **overnight** at 37°C with vigorous shaking.
3. In an Erlenmeyer culture flask, inoculate in 50 mL of LB medium + Resistance to the 1% (add 500 μ L of saturated culture) and incubate at 37 °C with vigorous shaking until the OD_{600 nm} is between 0.4 and 0.6.

Optimal transformation efficiencies are achieved when *E. coli* cells are at an optical density, OD_{600nm}, of between 0.4 and 0.6 (which will lead to 200 μ L of competent cells, using 50 μ L for each transformation). To obtain this biomass, the cell suspension should be **checked every 30 minutes** using a spectrophotometer. It usually takes between 2 and 3 hours after inoculation of the culture to reach this point.

In Sterility:

4. Move the content from the culture flask to a 50 mL Falcon tube and keep it on ice for 5-10 minutes.
5. Pellet the cells by centrifugation at 4600g for 7 minutes at **4°C** and discard the supernatant.
6. Resuspend the cells in 1 mL of ice-cold sterile distilled water. Add now 30 mL of ice-cold sterile distilled water and centrifugate again. It is important to **keep the cells on ice from this point forward** for a high transformation efficiency to be achieved. An increase in cell temperature will result in lower transformation efficiencies.

7. Repeat step 6 until the cells have been washed in water at least three times.

The cells must be washed extensively in sterile distilled water to remove the growth medium, which may not be suitable for electroporation; for example, it may have a high salt concentration, which will result in sample arcing.

8. Resuspend the cells in 1 mL of ice-cold 10% glycerol and move it to an 1,5 mL Eppendorf. Centrifuge at 10000g for 30 seconds at 4 °C.
9. Discard the supernatant and repeat step 8.
10. Resuspend the cells in a final volume of 200 µL of ice-cold 10% glycerol.
11. If cells are to be used immediately, place the cells on ice. To store electrocompetent cells, aliquot into 50 µL amounts and immediately snap-freeze using a liquid nitrogen bath and store at -80°C.