






## pET-28a(+) Backbone Amplification




### Introduction

This protocol allows to amplify the pET-28a(+) backbone using primers with Bsal restriction enzyme sites designed in Geneious.

### Reagents

-  pET-28a(+) backbone
-  pET-28a(+) backbone forward primer
-  pET-28a(+) backbone reverse primer
-  Q5 Hi-Fidelity 2X Master Mix
-  Milli-Q H<sub>2</sub>O

### Equipment

-  PCR tubes
-  Thermocycler
-  Pipette and tips

### Procedure

1. Add the following amounts of reagents to a sterile PCR tube:

Reagent	Amount
pET-28a(+) backbone (10 ng/μl)	1 μL
pET-28a(+) forward primer (10 μM)	1.25 μL
pET-28a(+) reverse primer (10 μM)	1.25 μL
Q5 Hi-Fi 2X master mix	12.5 μL
Milli-Q H <sub>2</sub> O	9 μL

2. Use the following PCR thermocycler program:

Step	Temperature	Time
Initial Denaturation	98°C	30 seconds
<b>Repeat the following 3 steps for 27 cycles</b>		
Denaturation	98°C	7 seconds
Annealing	72°C	15 seconds
Elongation	72°C	2 min 50 sec
Final Elongation	72°C	2 minutes
Hold	4°C	∞